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
Instruments for Treating
Condylar Process Fractures
acc. to Eckelt and Rasse

Fallbild von Dr. Hamid-Reza Sarajian, Rotes Kreuz Krankenhaus, Kassel, Deutschland
Case image provided by courtesy of Dr. Hamid-Reza Sarajian, Red Cross Hospital (RKK) Kassel, Germany

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Instruments for Treating Condylar Process Fractures acc. to Eckelt and Rasse

Introduction

The surgical treatment of mandibular condyle process fractures has seen significant advances in the past 20 years, thanks to important new contributions in the form of osteosynthesis techniques guaranteeing functional stability, such as the miniplate and the lag screw. Depending on fracture localization – diacapitular fractures, condylar neck fractures or basal condylar process fractures – special surgical techniques have been developed that differ both in terms of joint access and osteosynthesis technique used. Given a surgical indication, the surgeon therefore can select from a whole range of different techniques, each one requiring a special set of instruments for their safe application.


The selection of instruments presented here is intended to provide the surgeon with just the right instruments for the surgical techniques used.

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
Instruments

- ①+② Rasse hooklet
- ③ Tessier hooklet
 Microhooklets for repositioning and fixing dislocated fragments
- ④ Transbuccal set, modified acc. to Lauer/Feller, sharpened on one side, concave plate design for improved intraoral overview during screwdriver insertion

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


Surgical technique

The trocar is introduced through the cheek, then the trocar sleeve is pushed over the trocar through the cheek from inside the oral cavity. Following trocar removal, the retaining ring is applied to the trocar sleeve from the outside and the cheek compressed as a result. The concave form of the inner surface of the trocar sleeve plate gives the surgeon a much better view of the surgical site. The pictures illustrate screw insertion through the trocar.

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
Instruments

- 5 Set acc. to Rasse, with trocar and soft-tissue sleeve and matching circular notches for the grasping forceps no. 25-473-22-07 or 25-478-22-07.
- 6 Soft-tissue retraction forceps acc. to Rasse for the trocar system no. 25-473-00-07 or 25-478-00-07, for intra-oral application.
- 7 Soft-tissue retraction forceps acc. to Rasse for the transbuccal set no. 25-478-00-07, for use with preauricular access.
- 8 Transbuccal set for resorbable systems acc. to Rasse, with matching circular notches for the grasping forceps no. 25-479-22-07.
- 9 Soft-tissue retraction forceps for the transbuccal set no. 25-479-00-07, for preauricular access.

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

Surgical technique

This transbuccal set represents an alternative to the surgical set no. 25-473-00-07.

The soft-tissues are held off with the Rasse soft-tissue retraction forceps applied to the trocar from oral.

When using preauricular access, it may be necessary in certain cases to use a transbuccal trocar for lag screw insertion, depending on the course of the fracture. The soft-tissue retraction forceps applied to the trocar system from the inside ensures a good view of the surgical site during the intervention.

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Ø 4 mm (inner diameter)
for resorbable screws

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Instruments

- ⑩ + ⑪ Periosteal elevators, for dissection and elevation in
very deep regions of the mandibular joint, especially
with submandibular/perimandibular access.

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
Instruments

- ⑫ + ⑬ Periosteal elevator, extra long, slightly and strongly curved. For dissection and elevation in very deep regions at the rear edge of the mandibular ramus when using oral access.

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
Instruments

- 14 Eckelt retractor
 Retractor with wide-open hook for easier mandibular application when using oral access, modified acc. to Eckelt.
- 15 Rasse retractor
 Retractor with wide-open hook for easier mandibular application when using oral access, modified acc. to Rasse.

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
Instruments

- 16 Eckelt retractor with cold light
Retractor with wide-open hook for easier mandibular application when using oral access, with cold light for surgical site illumination.

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
Instruments

- 17 Schuchardt/Lauer progeny retractor
External progeny retractor modified acc. to Lauer, with window, for use with oral access and trans-buccal screw insertion.
- 18 Metz/Lauer progeny retractor
Internal progeny retractor modified acc. to Lauer, with window, for use with oral access and trans-buccal screw insertion.

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
Surgical technique

When using oral access, the window provided in the Schuchardt/Lauer and Metz/Lauer retractor allows the surgeon to insert screws for miniplate fixation while simultaneously retracting the cheek soft-tissues.

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
Instruments

- 19 Schuchardt/Lauer progeny retractor with cold light
External progeny retractor modified acc. to Lauer, with window, for use with oral access and transbuccal screw insertion.
With additional cold-light connector for better illumination of the surgical site.

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
Instruments

- 20 Metz/Lauer progeny retractor with cold light
Internal progeny retractor modified acc. to Lauer, with window, for use with oral access and transbuccal screw insertion.
With additional cold-light connector for better illumination of the surgical site.

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
Instrumente

- 21 Schuchardt/Lauer progeny retractor
External progeny retractor modified acc. to Lauer, with window, for use with oral access and transbuccal screw insertion.
With two additional guide tubes for fitting an optical system.

- 22 Metz/Lauer progeny retractor
Internal progeny retractor modified acc. to Lauer, with window, for use with oral access and transbuccal screw insertion.
With two additional guide tubes for fitting an optical system.

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
Surgical technique

When using oral access, the window provided in the Schuchardt/Lauer retractor allows the surgeon to insert screws for miniplate fixation while simultaneously retracting the cheek soft-tissues.

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

Instruments

- ②③ + ②④ Long instrument for repositioning and manipulating dislocated and fractured mandibular joints.
- ②⑤ Eckelt repositioning pin
 Pin for repositioning dislocated fragments via a transbuccal trocar system (see ④ and ⑤).
 Requires handle no. 25-402-99-07. Can be used irrespective of the type of surgical access.

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
Surgical technique

The repositioning pin can be screwed into the small fragment via the transbuccal trocar system (25-478-00-07) with any surgical access – preauricular, periangular/submandibular or intraoral. This facilitates bone repositioning and fixation during the osteosynthesis procedure.

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
Instruments

- 27 Eckelt plate-pulling instrument
Plate-holding and plate-pulling instrument for intraoral support of the osteosynthesis plate during screw insertion, as well as for pulling the plate at its distal end after fitting it to the proximal fragment, thus moving it caudally to ensure precise adaptation of the fragments. Only for caudal pull.

Surgical technique

To enable exact adaptation of the fragments when using oral access, the miniplate is first fixed to the proximal fragment with two screws, then the plate-pulling instrument is inserted into the caudal-most free hole of the miniplate. Now the plate can be drawn caudally with the instrument for precise adaptation of the fragments and insertion of the third screw next to the fracture cleft can be effected.

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
Instruments


28 + 29 Eckelt repositioning forceps (left and right versions)

Modified repositioning forceps with particularly slender
jaws and optimized tooth arrangement. For repositioning
the proximal fragment when using extraoral
(perimandibular/submandibular) access.

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
Instruments

- 30 + 31 Eckelt/Rasse fracture forceps
 For moving the mandible caudally and supporting the repositioning process.


Surgical technique

Applied transcutaneously, the forceps are fixed to the mandibular angle. They are available in a right and left version (for right/left side of the mandible). The jaws or tongs of the forceps have a strongly curved (concave) shape to keep the cheek soft tissues out of the way.

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


All instruments at a glance


No. of Ill.	Item No.	Designation
1	15-006-10-07	Bone hooklet acc. to Rasse
2	15-068-05-07	Bone hooklet sharp acc. to Rasse
3	48-250-19-07	Skin hooklet sharp acc. to Tessier
4	25-473-00-01	Transbuccal set modified acc. to Lauer/Feller
5	25-478-00-07	Transbuccal set acc. to Rasse
	25-478-15-07	Trocar only for 25-478-00-07
	25-478-20-07	Drill sleeve only for 25-478-00-07
6	25-473-22-07	Soft tissue retraction forceps for 25-473-00-07
7	25-478-22-07	Soft tissue retraction forceps for 25-478-00-07
8	25-479-00-07	Transbuccal set acc. to Rasse
	25-479-15-07	Trocar only for 25-479-00-07
	25-479-20-07	Drill sleeve only for 25-479-00-07
9	25-479-22-07	Soft tissue retraction forceps for 25-479-00-07
10	38-633-01-07	Periosteal elevator slightly curved acc. to Krenkel
11	38-633-02-07	Periosteal elevator strongly curved acc. to Krenkel
12	37-544-01-07	Periosteal elevator slightly curved acc. to Krenkel
13	37-544-02-07	Periosteal elevator strongly curved acc. to Krenkel
14	38-502-24-01	Retractor acc. to Eckelt
15	48-121-20-07	Retractor acc. to Rasse
16	38-504-27-01	Retractor with cold light acc. to Eckelt
17	38-616-10-07	External progeny retractor acc. to Schuchardt/Lauer
18	38-616-20-07	Internal progeny retractor acc. to Metz/Lauer
19	38-618-10-07	External progeny retractor with cold light acc. to Schuchardt/Lauer
20	38-618-20-07	Internal progeny retractor with cold light acc. to Metz/Lauer
21	38-620-10-07	External progeny retractor with att. for optics acc. to Schuchardt/Lauer
22	38-620-20-07	Internal progeny retractor with att. for optics acc. to Metz/Lauer
23	38-633-03-07	Condylar repositioning instrument Fig. 11
24	38-633-04-07	Condylar repositioning instrument Fig. 12
25	25-430-20-07	Repositioning pin for handle 25-402-99-07 acc. to Eckelt
26	25-402-99-07	Standard handle only for 25-430-20-07
27	38-715-21-07	Plate-pulling instrument acc. to Eckelt
28	50-241-11-07	Repositioning forceps right acc. to Eckelt
29	50-241-12-07	Repositioning forceps left acc. to Eckelt
30	50-325-11-07	Fracture forceps left acc. to Eckelt/Rasse
31	50-325-12-07	Fracture forceps right acc. to Eckelt/Rasse

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Endoscopy-assisted open treatment of condylar fractures
of the mandible: extraoral vs. intraoral approach
Int J Oral Maxillofac Surg 2002; 31:237-243

Undt G, Kermer C, Rasse M, Sinko K, Ewers R:

Transoral miniplate osteosynthesis of condylar neck fractures
Oral Surg Oral Med Oral Pathol Oral Radiol Endod 1999;
88(5):534-543

KLS Martin Group

KLS Martin Australia Pty Ltd.

Sydney · Australia
Tel. +61 2 9439 5316
australia@klsmartin.com

KLS Martin do Brasil Ltda.

São Paulo · Brazil
Tel. +55 11 3554 2299
brazil@klsmartin.com

KLS Martin Medical (Shanghai) International Trading Co., Ltd.

Shanghai · China
Tel. +86 21 5820 6251
china@klsmartin.com

KLS Martin India Pvt Ltd.

Chennai · India
Tel. +91 44 66 442 300
india@klsmartin.com

Martin Italia S.r.l.

Milan · Italy
Tel. +39 039 605 67 31
italia@klsmartin.com

Nippon Martin K.K.

Tokyo · Japan
Tel. +81 3 3814 1431
nippon@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Penang · Malaysia
Tel. +604 505 7838
malaysia@klsmartin.com

Martin Nederland/Marned B.V.

Huizen · The Netherlands
Tel. +31 35 523 45 38
nederland@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Moscow · Russia
Tel. +7 499 792-76-19
russia@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

KLS Martin UK Ltd.

London · United Kingdom
Tel. +44 1189 000 570
uk@klsmartin.com

KLS Martin LP

Jacksonville · Florida, USA
Tel. +1 904 641 77 46
usa@klsmartin.com

Gebrüder Martin GmbH & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · 78532 Tuttlingen · Germany
P.O. Box 60 · 78501 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com