






BRM[®]
Extremities

Digitalis

5 Metacarpophalangeal implants
4 Proximal interphalangeal implants



MCP



PIP


Osteoarthritis
Rheumatoid arthritis
Destroyed articular surfaces
Ankylosed joints or with limited range of motion
Non-functional joint due to inadequate bony alignment



IMPLANT

CODE

DESCRIPTION



| | |
|------------|-------------------------------|
| DDG3T02001 | DIGITALIS PIP Spacer – Size 1 |
| DDG3T02002 | DIGITALIS PIP Spacer – Size 2 |
| DDG3T02003 | DIGITALIS PIP Spacer – Size 3 |
| DDG3T02004 | DIGITALIS PIP Spacer – Size 4 |

| | |
|------------|-------------------------------|
| DDG3T01001 | DIGITALIS MCP Spacer – Size 1 |
| DDG3T01002 | DIGITALIS MCP Spacer – Size 2 |
| DDG3T01003 | DIGITALIS MCP Spacer – Size 3 |
| DDG3T01004 | DIGITALIS MCP Spacer – Size 4 |
| DDG3T01005 | DIGITALIS MCP Spacer – Size 5 |


Main contraindications

- Inadequate musculo-tendon and skin system
- Inadequate neuro-vascular system
- Bone demineralization at a significant stage
- Inadequate bone stock
- Child patient
- Infection in progress and active sepsis

INSTRUMENTS

CODE

DESCRIPTION



| | |
|-----------------|---|
| SET.B.DIGITALIS | Digitalis Instruments Set without trials |
| DDGI000000B | Digitalis Small sterilization box (empty) |
| DDGI203000 | Digitalis Handle |
| DDGI201003 | Digitalis MCP Multisize Trial |
| DDGI201004 | Digitalis PIP Multisize Trial |
| DDGI201001 | Digitalis MCP Reamer |
| DDGI202001 | Digitalis MCP Rasp |
| DDGI201002 | Digitalis PIP Reamer |
| DDGI202002 | Digitalis PIP Rasp |
| SET.DIGITALIS | Digitalis complete Instruments Set |
| DDGI000001B | Digitalis instruments sterilization box (empty) |
| DDGI201013 | Silicone MCP Digitalis trial - size 1 |
| DDGI201023 | Silicone MCP Digitalis trial - size 2 |
| DDGI201033 | Silicone MCP Digitalis trial - size 3 |
| DDGI201043 | Silicone MCP Digitalis trial - size 4 |
| DDGI201053 | Silicone MCP Digitalis trial - size 5 |
| DDGI201014 | Silicone PIP Digitalis trial - size 1 |
| DDGI201024 | Silicone PIP Digitalis trial - size 2 |
| DDGI201034 | Silicone PIP Digitalis trial - size 3 |
| DDGI201044 | Silicone PIP Digitalis trial - size 4 |